

Weekly Eyewash Inspection Log

Location_____

Year: 2026

Instructions: Activate emergency eyewash for at least 1 minute to ensure proper operation:

- Ensure the dust caps pop off when eyewash is activated
- Ensure streams of water meet in the middle/sufficiently can reach user eyes of the eyewash.
- Ensure the water flows in 1 second or less of activation and remain on until it is intentionally turned off
- Ensure water runs clear and is visibly free from foreign particles
- Ensure access to emergency eyewash is unobstructed. Ensure eyewash is clean and dust caps are in place to prevent dust from entering the eyewash.
- If there are any issues, tag unit out of service and enter a facility work order to request that the unit be repaired as soon as possible. Keep weekly eyewash inspection log sheet for 3 years.

Week of	Date	Initials	Week of	Date	Initials	Week of	Date	Initials
12/28/2025			5/3/2026			9/6/2026		
1/4/2026			5/10/2026			9/13/2026		
1/11/2026			5/17/2026			9/20/2026		
1/18/2026			5/24/2026			9/27/2026		
1/25/2026			5/31/2026			10/4/2026		
2/1/2026			6/7/2026			10/11/2026		
2/8/2026			6/14/2026			10/18/2026		
2/15/2026			6/21/2026			10/25/2026		
2/22/2026			6/28/2026			11/1/2026		
3/1/2026			7/5/2026			11/8/2026		
3/8/2026			7/12/2026			11/15/2026		
3/15/2026			7/19/2026			11/22/2026		
3/22/2026			7/26/2026			11/29/2026		
3/29/2026			8/2/2026			12/6/2026		
4/5/2026			8/9/2026			12/13/2026		
4/12/2026			8/16/2026			12/20/2026		
4/19/2026			8/23/2026			12/27/2026		
4/26/2026			8/30/2026					