

## MANAGEMENT & DISPOSAL METHODS FOR YALE UNIVERSITY CLINIC CHEMICAL WASTES

<u>CHEMICAL/DRUG/PRODUCT NAME</u>	<u>MANAGEMENT/DISPOSAL</u>	<u>HAZARD CLASS</u>
3M Attest - G. Stearothermophilus	APW	
3M No Sting Liquid Bandage	HW	Flammable
3M Remover Lotion	HW	Flammable
5-Fluorouracil	HW	
Abacavir	APW	
Acamprosate Calcium	APW	
Acellular Pertussis Vaccine	APW	
Acetaminophen	APW	
Acetasol (Acetic Acid & Hydrocortisone)	APW	
Acetic Acid 1-5% / Water Solutions	D	
Acetic Acid 2% in Aqueous Aluminum Acetate	N-HW	
Acetone	HW	Flammable
Activate 10:1 Bleach	HW	Corrosive (Base)
ACU-DTM Dermatophyte Test Medium	APW	
AddiPaK	APW	
Adenosine	APW	
Adhesive Tape Remover Pads (UNUSED)	HW	Flammable
Ajax Oxygen Bleach Cleaner	N-HW	
AK-DILATE Phenylephrine Ophthalmic Solution	APW	

Albuterol	APW	
Alcare-Antiseptic Foam	HW	Flammable
Alcohol Prep Pads/Swabs (UNUSED)	HW	Flammable
Alendronate	APW	
Aliskiren	APW	
Almotriptan Malate	APW	
Aluminum Chloride Crystals	N-HW	
Aluminum Chloride/Water Solution	N-HW	
Alteplase	APW	
Amdoxovir	APW	
Amerse 2	HW	Corrosive (Base)
Amlodipine	APW	
Amlodipine Besylate	APW	
Ammonia Inhalant (UNUSED)	HW	Flammable
Ammonia Spirits	HW	Flammable
Amoxicillin	APW	
Amphotericin B	APW	
Ampicillin	APW	
Amprenavir	APW	
Aprepitat	APW	
Apricitabine	APW	
Aptima Unisex Swab Specimen Collection Kit (unused)	N-HW	
Aqua Shine Aerosol	HW	Flammable
Aqua Treat Aerosol	HW	Flammable
Aquasonic 100 Ultrasound Gel	T	
Aripiprazole	APW	
Arsenic Trioxide & its packaging	HW	Toxic
Ascend Germicidal Detergent	N-HW	
Asenapine Maleate	APW	

Asparaginase Powder	APW	
Aspirin	APW	
Atazanavir	APW	
Atomoxetine	APW	
Atorvastatin Calcium	APW	
Atovaquon	APW	
Atropine	APW	
Atropine Sulfate	APW	
Avagard D	HW	Flammable
Azaserine	HW	Toxic
Azorvastatin Calcium	APW	
Bacitracin Ointment	APW	
Balsalazide	APW	
Barium Sulfate Suspension	HW	Toxic
Baserunner	HW	Corrosive (Base)
BD Bactec Culture Vials	APW	
Becaplermin	APW	
Benazepri	APW	
Benzalkonium Chloride Towelettes (UNUSED)	APW	
Benzoin Compound Tincture	HW	Flammable
Benzoin Crystals	N-HW	
Benzoin Tincture	HW	Flammable
Bevacizumab	APW	
BiCNU (vial with ethanol)	HW	Flammable
Biozide	HW	Flammable
Birth Control Pills	APW	
Bleach	HW	Corrosive (Base)
Bleomycin Sulfate	APW	
Boric Acid Powder	T	

Botox	APW	
Botulinum Toxins	APW	
Budesonide	APW	
Buprenorphine	Controlled Substance	
Bupivacaine	APW	
Bupropion	APW	
Buspiron	APW	
Butyl Cellosolve	HW	Flammable
Butyric Acid including n-Butyric Acid	N-HW	
Caffeine Citrate	APW	
Candela Gentle Cool Cryogen	N-HW	
Candesartan Cilexetil	APW	
Cantharone Plus Wart Treatment	HW	Flammable/Corrosive (Acid)
Carbamazepine	APW	
Carbasone	HW	Toxic
Carbidopa	APW	
Carbolic Acid/Phenol	HW	Toxic
Carboplatin	APW	
Carboxymethyl Cellulose Powder	N-HW	
Castile Soap Towelettes (UNUSED)	T	
Cavicide	N-HW	
Caviwipes (UNUSED)	N-HW	
Cefolinir	APW	
Cefotaxime	APW	
Cefprozil	APW	
Ceftrazidime	APW	
Ceftriaxone	APW	
Cell-Dyn 3200 WBC Lyse	D	
Centrorelix Acetate	APW	

Cephalexin	APW	
Cerefolin	APW	
Cetacaine Spray	HW	Flammable
Cevimeline	APW	
Chemstrips	T	
Chemtronics Freeze Spray	N-HW	
Chlorambucil	HW	Toxic
Chloraprep	HW	Flammable
Chlorazol Black Powder	N-HW	
Chlorazol Black E - Delasco	HW	Corrosive (Base)
Cholesterol Powder	N-HW	
Choriogonadotropin Alfa	APW	
Cidex Plus 28 Day Solution	N-HW	
Cidex Activated Dialdehyde Solution	D	
Cidex OPA Solution	D	
Cilostazol	APW	
Ciprofloxacin	APW	
Cisplatin	APW	
Citalopram	APW	
Citrace	HW	Flammable
Citrus II CPAP Mask Cleaner	D	
Cladribine	APW	
Clindamycin – except Topical Solution	APW	
Clindamycin Topical Solution	HW	Flammable
Clopidogrel	APW	
Clorox Disinfecting Spray	HW	Flammable
Clorox Disinfecting Wipes (UNUSED)	N-HW	
Coal Tar Topical Solution	HW	Flammable
Coamatic (Substrate S02403 vials only)	HW	Toxic

Cobalamin	APW	
Coil Disinfectant	N-HW	
Colesevelam	APW	
Collagen	APW	
Collodion	HW	Flammable
Collodion Flexible	HW	Flammable
Collodion Remover	HW	Flammable
Comet Cleaner w/ Bleach - Liquid	HW	Corrosive (Base)
Comet Cleaner w/ Bleach - Powder	N-HW	
Comet Disinfecting Bathroom Cleaner	N-HW	
Comet Powder	N-HW	
Coumadin tablets & their packaging	HW	Toxic
Crew Bathroom Cleaner and Scale Remover #44	N-HW	
Cry-Ac-3 Liquid Nitrogen	N-HW	
Cutex Non-Acetone Nail Polish Remover	HW	Flammable
Cyclophosphamide	HW	Toxic
Cytarabine	APW	
Cyto Fix	HW	Flammable
CytoRich Red	APW	
CytoRich Collection Fluid	HW	Flammable
D(+) Xylose	T	
D5W Dextrose/Water	D	
Dacarbazine	APW	
Dactinomycin	APW	
Dalteparin	APW	
DAP Kwik Foam Polyurethane Foam	HW	Flammable
Darbepoetin Alfa	APW	
Darifenacin	APW	
Darunavir	APW	

Daunorubicin	HW	Toxic
Deep Sea Nose Spray	APW	
DeFogger	N-HW	
Delavirdine Mesylate	APW	
Dermabond (UNUSED)	N-HW	
Dermabond (USED)	APW	
Derma-Smoothe/FS	N-HW	
Desloratadine	APW	
Desmopressin Acetate	APW	
Dexamethasone	APW	
Dexrazoxane	APW	
Dextrose	D	
Dextrose/Sodium Chloride Solution	D	
Diazepam	Controlled Substance	
Digestive Advantage Medical Food	APW	
Dihydroergotamine Mesylat	APW	
Diltiazem	APW	
Dimethyl Acetylenedicarboxylate	HW	Corrosive (Acid)
Dimethyl Sulfoxide/Water/Lidocaine/Dexamethosone Sodium Phosphate/Heparin	D	
Diphtheria + Tetanus Toxoids	APW	
Dishwashing Soaps - Liquid	D	
Dispatch Hospital Cleaner Disinfectant with Bleach	N-HW	
Divalproex Sodium	APW	
Dolasetron	APW	
Donepezil	APW	
Doxorubicin	APW	
Doxycycline	APW	
Drospirenone	APW	
Drysol - Person & Covey	HW	Flammable

Dulbelcco's Medium	D	
Duloxetine	APW	
Dust-Off Cleaning Wipes	T	
DynaPine	HW	Flammable
Dynaquat Cleaner	N-HW	
Ecolab Gold Label Mist	HW	Flammable
Efavirenz	APW	
Elite Professional Bleach	HW	Corrosive (Base)
Elmer Multipurpose Glue	T	
Elvitegravir	APW	
Elvucitabine	APW	
Emtricitabine	APW	
Endozime AW Plus	N-HW	
Endozime SLR	N-HW	
Enfuvirtide	APW	
Enoxaparin Sodium	APW	
Entacapone	APW	
Envirocide	N-HW	
Ephedrine Sulfate	APW	
Epinephrine & its packaging	HW	Toxic
Epinephrine/Bupivacaine mixtures	APW	
Epinephrine/Lidocaine mixtures	APW	
Epinephrine/Marcaine mixtures	APW	
Epinephrine/Xylocaine mixtures	APW	
Epoetin Alfa	APW	
Erythromycin - <b>except 2% Gel</b>	APW	
Erythromycin Gel, 2%	HW	Flammable
Escitalopram	APW	
Esomeprazole Magnesium	APW	

Estradiol	APW	
Estrogen	APW	
Ethanol Solution	HW	Flammable
Ethinyl Estradiol	APW	
Etoposide	HW	Flammable
Etravirine	APW	
Ezetimibe	APW	
Expo Cleaner for Dry Erase White Boards	HW	Flammable
Exsel Shampoo	HW	Toxic
Eyewash Solution/Irrigation	D	
Fabulous Cleaner	HW	Flammable
Fade-A-Dyne II	HW	Flammable
Famotidine	APW	
Fantastik Cleaner	D	
Fenofibrate	APW	
Ferric Chloride/Water	HW	Corrosive (Acid)
Ferric Subsulfate Solution - Monsel's Solution	HW	Corrosive (Acid)
Fisher Scientific Acid Spill Cleanup	N-HW	
FIST Disinfecting Bowl Cleaner	HW	Corrosive (Acid)
Fleet Enema	D	
Fluconazole	APW	
Flumazenil	APW	
Fluorescein	N-HW	
Fluorescein Injection	N-HW	
Fluoxetine	APW	
Fluphenazine Decanoate	APW	
Fluticasone Propionate	APW	
Follitropin Alfa	APW	
Follitropin Beta	APW	

Formaldehyde, 37% Solution	HW	Flammable/Toxic
Formalin (UNUSED)	HW	Toxic
Formalin (USED)	N-HW	
Formula 409 All Purpose Cleaner	N-HW	
Formula 409 Glass & Surface Cleaner	N-HW	
Fosamprenavir Calcium	APW	
Frovatriptan Succinate	APW	
Galantamine	APW	
Ganciclovir	APW	
Ganirelix Acetate	APW	
Gas-X	APW	
GBG Aloe Gel Instant Hand Sanitizer w/ Aloe	HW	Flammable
Granulocyte-Colony Stimulating Factor	APW	
Gebauer's Fluori-Methane	N-HW	
Gebauer's Pain Ease	N-HW	
Gebauer's Ethyl Chloride	HW	Flammable
Generation Block	N-HW	
Genetian Violet	APW	
Genuine EC2 Electrode Paste	T	
Genzyme OSOM Card Pregnancy Test	N-HW	
Glade Aerosol Sprays	HW	Flammable
Glance NA Cleaning Product	N-HW	
Glucagon	APW	
Glucola	D	
Glucose Tablets	APW	
Glutaraldehyde 2.6%/Water Solution	D	
Glute Out	N-HW	
Glute Out -2	N-HW	
Glut-Rx	N-HW	

Glycerin USP	N-HW	
Glycolic Acid/Water Solution - Delasco	HW	Corrosive (Acid)
Gold Label Aerosol Spray	HW	Flammable
Golytely Disposable Jug	N-HW	
Goof Off 2000	HW	Flammable
Googone	HW	Flammable
GP Forward SC	N-HW	
Granisetron	APW	
Haemo-Sol	T	
Hair Gels	D	
Haloperidol	APW	
Hand Soaps/Cleansers - Liquid	D	
hCG Pregnancy Test Strips	T	
Hemoccult Developer	HW	Flammable
Heparin Sodium Injection	APW	
Hepatitis Vaccines	APW	
Hibiclens	D	
Histamine Phosphate/Glycine	APW	
Histofreezer	HW	Flammable
Hospicide	HW	Flammable
Hospiseptic Surface Disinfectant	HW	Flammable
Hot Shot Spider Killer	HW	Flammable
Human Chorionic Gonadotropin (hCG)	APW	
Hurricane Topical Anesthetic Spray	HW	Flammable
Hurricane Spray	HW	Flammable
Hycamtin	APW	
Hydrochloric Acid Solution	HW	Corrosive (Acid)
Hydrochlorothiazide	APW	
Hydrocortisone Cream/Ointment	APW	

Hydrogel Wound Dressing	T	
Hydrogen Peroxide 3%/Water Solution	D	
Hydroxyurea	APW	
Hydroxyzine	APW	
Hyoscyamine	APW	
Ibuprofen	APW	
IC-GREEN Indocyanine Green Dye (Powder & in Aqueous Solvent)	N-HW	
IC-GREEN Aqueous Solvent	D	
Idarubicin	HW	Toxic
Iloperidone Tablets	APW	
Immune Globulin	APW	
Immunitin	APW	
Indermil Tissue Adhesive	N-HW	
Indinavir Sulfate	APW	
Infliximab	APW	
Influenza Vaccine	APW	
Instant Hot Compresses	T	
Instat Collagen Absorbable Hemostat	APW	
Insulin preserved with m-Cresol	HW	Toxic
Insulin without m-Cresol	APW	
Interferon	APW	
Iodine Pads/Swabsticks (UNUSED)	T	
Iodine Tincture	N-HW	
Iodixanol	N-HW	
Iodoform Packing Strips	N-HW	
Irbesartan	APW	
Irbesartan - Hydrochlorothiazide	APW	
Isoniazid	APW	
Isopropanol/Isopropyl Alcohol, 70% Solution	HW	Flammable

Isopto Homatropine	APW	
Isradipine	APW	
Kimwipes Lens Cleaning Solution	N-HW	
Klenzyme	N-HW	
Kodak Intensifying Screen Cleaner & Antistatic Solution	N-HW	
Kwell Shampoo	HW	
KY Jelly	T	
Lacrilube	APW	
Lamivudine	APW	
Lanzoprazole	APW	
Latanoprost Ophthalmic Solution	APW	
Leuprolide Acetate	N-HW	
Levetiracetam	APW	
Levodopa	APW	
Levofloxac	APW	
Levomilnacipran	APW	
Levulan Kerasticks	APW	
Lidocaine (syringe, cream, ointment, patch, etc...)	APW	
Lidocaine + Epinephrine (syringe, cream, ointment, etc...)	APW	
Lidocaine + Prilocaine (cream)	APW	
Linaire Gel	N-HW	
Lindane	HW	Toxic
Linen Spray	HW	Flammable
Liquasonic Ultrasound Gel	N-HW	
LiquiBand (UNUSED)	N-HW	
LiquiBand (USED)	APW	
Liquid Paper White Out	HW	Flammable
Lissamine Green Ophthalmic Strips	APW	
Loratadine	APW	

Lorazepam	Controlled Substance	
Losartan Potassium	APW	
Lovastatin	APW	
Lubricating Grease	N-HW	
Lugol's Solution	N-HW	
Lurasidone Hydrochloride	APW	
Lysol All Purpose Cleaner 4 in 1	N-HW	
Lysol Disinfectant Spray	HW	Flammable
Lysol spray	HW	Flammable
Maalox	D	
Mannitol	APW	
Maraviroc	APW	
Mastisol Liquid Adhesive	HW	Flammable
Mavidon Medical Adhesive Remover	HW	Flammable
Maxforce FC Roach & Ant Baits	N-HW	
Maxipime	APW	
McKesson Instant Hand Sanitizer	HW	Flammable
McKesson Lubricating Jelly	T	
Mechloroethamine	APW	
MediChoice Lubricating Jelly	T	
MediChoice Ultrasound Gel	T	
Medi-Pak Instant Hand Sanitizer	HW	Flammable
Medisense Glucose & Ketone Control Solution	N-HW	
Medline Packing Strips	APW	
Medovations WeightRight Bougies (UNUSED)	N-HW	
Medphene Plus Spray	HW	Flammable
Melphalan	HW	Toxic
Memantine	APW	
Meningococcal Vaccine	APW	

Menotropins for injection	APW	
Menthol Crystals	N-HW	
Mepivacaine	APW	
Mercury Containing Devices (i.e. Thermometers, Gauges, etc...)	HW	Toxic
Mercury, Ammoniated	HW	Toxic
Mesalamine	APW	
Methenamine	APW	
Methotrexate	APW	
Methylene Blue Crystals	N-HW	
Methylene Blue/Water	D	
Metoprolol	APW	
Metricide - Activated	D	
Metricide - Unactivated	N-HW	
Metrizyme	N-HW	
Metronidazole	APW	
Michel's Fixative	N-HW	
Microban Disinfectant Spray	N-HW	
Microban QGC	N-HW	
Midazolam	Controlled Substance	
Mineral Oil	N-HW	
Mirtazapine	APW	
Mitomycin C	HW	Toxic
Mitoxantrone	APW	
Modafinil	Controlled Substance	
Monsel's-Ferric Subsulfate Solution (EMPTY CONTAINERS ONLY)	APW	
Monsel's-Ferric Subsulfate Solution (NON-EMPTY CONTAINERS)	HW	Corrosive (Acid)
Montelukast	APW	
Morphine	Controlled Substance	
Motrin	APW	

Moxifloxacin	APW	
Mupirocin	APW	
Murphy's Oil Soap	D	
Nafcillin	APW	
Naloxone	APW	
Naproxen	APW	
Nebivolol	APW	
Nefazodone	APW	
Nelfinavir Mesylate	APW	
Niacin	APW	
Nicotine & its packaging	HW	Toxic
Nisoldipine	APW	
Nitroglycerine & its packaging	HW	Toxic
Norgestimate	APW	
Normal Saline - in bottles & IV bags	D	
Normal Saline - in syringes	APW	
Novus Fine Scratch Remover #2	N-HW	
Nu-Gauze	APW	
NuPrep Skin Prep Gel	T	
Nystatin	N-HW	
Odor Eliminator	N-HW	
Olanzapine	APW	
Olmesartan Medoxomil	APW	
Omega 3 Acid Ethyl Esters	APW	
Ominpaque Injection	D	
Omni Prep Paste	T	
Ondansetron	APW	
Orthozime	HW	Flammable
Oust Air Sanitizer	HW	Flammable

Oust Disinfectant Spray	HW	Flammable
Oxacillin	APW	
Oxcarbazepine	APW	
Oxiconazole Nitrate Cream	APW	
Oxybutynin Chloride	APW	
Paclitaxel Injection	HW	Flammable
Paliperidone Palmitate	APW	
Palonosetron	APW	
Paroxetine	APW	
PAWS Wipes (UNUSED)	HW	Flammable
PDI Sani-Coth Plus Wipes (UNUSED)	N-HW	
PDI Super Sani-Cloth Wipes (UNUSED)	HW	Flammable
Pegaspargase	APW	
Pegfilgrastim	APW	
Pentamidine	APW	
Peppermint Oil	N-HW	
Peppermint Spirit	HW	Flammable
Perindopril Erbumine	APW	
Phenol Powder & Solution	HW	Toxic
Phenol Swabs (UNUSED)	HW	Toxic
Phenol Swabs (USED)	APW	
Phenylmercuric Acetate preserved medications	HW	Toxic
Phisoderm	D	
PhisoHex	HW	Toxic
Phosphate Buffered Saline	D	
Physostigmine & its packaging	HW	Toxic
Physostigmine Salicylate & its packaging	HW	Toxic
Phytonadione	APW	
Pine Oil Cleaner	N-HW	

Pitney Bowes EZ Seal Solution	N-HW	
PMOP All Purpose Correction Fluid	HW	Flammable
Pneumococcal Vaccine	APW	
Podophyllum Resin	N-HW	
Podophyllum Resin in Tincture Benzoin - Delasco	HW	Flammable
Porcine Xenograft - <b>not in a chemical preservative</b>	APW	
POR-SAN	N-HW	
Potassium Bicarbonate Powder	T	
Potassium Hydroxide 10%/Water Solution	HW	Corrosive (Base)
Povidone Iodine Prep Pads & Swabsticks (UNUSED)	T	
Povidone Iodine Solution	N-HW	
PPD Test	APW	
Pramipexole Dihydrochloride	APW	
Pramoxine	APW	
Pravastatin Sodium	APW	
Prednisone	APW	
Pregabalin	Controlled Substance	
PreservCyt Solution	HW	Flammable
Progesterone	APW	
Progesterone Gel	APW	
PRO STRIP	HW	Corrosive (Base)
Proleukin	APW	
Proparacaine Ophthalmic Solution	APW	
Propylene Glycol Liquid	N-HW	
ProsALL Neutra Odor	HW	Flammable
Protopic Tacrolimus Ointment	APW	
PSS Select Disinfecting/Deodorizing/Cleaning Wipes (UNUSED)	HW	Flammable
Purell Hand Sanitizer	HW	Flammable
Pyrodoxine	APW	

Quanto Disinfectant	N-HW	
Quetiapine Fumarate	APW	
Quick Shine	HW	Flammable
Quick Vue Strep A Test	APW	
Quik-Dip Stain Solution 1	HW	Flammable/Toxic
Quik-Dip Stain Solution 2	N-HW	
Quik-Dip Stain Solution 3	N-HW	
Racivir	APW	
Raltegravir	APW	
Ramipiril	APW	
Ranitidine	APW	
Ranolazine	APW	
Reagent Alcohol	HW	Flammable
Red Z	N-HW	
Remove Wipes (UNUSED)	HW	Flammable
Remune	APW	
Renografin-60	N-HW	
Reserpine	HW	Toxic
Resorcinol	HW	Toxic
Rifaximin	APW	
Rilpivirine	APW	
RIMSO-50	N-HW	
Risperidone	APW	
Ritonavir	APW	
Rituximab	APW	
Rivastigmine Tartrate	APW	
Ropinirole	APW	
Rosuvastatin Calcium	APW	
Safetex Scientific Cytology Fixative	HW	Flammable

Salicylic Acid Powder	T	
Sani-Cloth Bleach Wipes (UNUSED)	HW	Corrosive (Base)
Sani-Cloth Disinfecting Wipes (UNUSED)	N-HW	
Sani-Cloth Plus Wipes (UNUSED)	N-HW	
Sani-Dex Alc Wipes (UNUSED)	HW	Flammable
Saquinavir Meslyate	APW	
Scrub Free Soap Scum Remover	HW	Corrosive (Acid)
Selsun Shampoo	HW	Toxic
Seracult Developer	HW	Flammable
Shampoos/Conditioners - except Kwell, Exsel, and Selsun	D	
SHINE UP Lemon Furniture Polish	HW	Flammable
Signa Creme	T	
Signa Spray	D	
Silver Nitrate Applicators/Sticks (UNUSED)	HW	Oxidizer/Toxic
Silver Nitrate Crystals	HW	Oxidizer/Toxic
Silver Nitrate Solution	HW	Oxidizer/Toxic
Silver Sulfadiazine Cream	HW	Toxic
Simoniz Disinfectant Cleaner Spray	HW	Flammable
Simoniz Glass Cleaner Spray	HW	Flammable
Simple Green	D	
Simply Clean	N-HW	
Simvastatin	APW	
Sinografin	N-HW	
Skin Moisturizers/Hand Creams	APW	
Sklar Instruments Ready to Use Paste	APW	
Sodium Bicarbonate 8.4%/Water Solution	D	
Sodium Bicarbonate Injection	APW	
Sodium Bicarbonate Powder	T	
Sodium Chloride - in syringes	APW	

Sodium Chloride Solution - in bottles & IV bags	D	
Sodium Meta-Bisulfite Granular	N-HW	
Sodium Nitrite - in Strep Kit	HW	Oxidizer
Sodium Oxybate	APW	
So-Drox	HW	Corrosive (Base)
Solifenacin Succinate	APW	
Sorbitol	N-HW	
Special Formula Shredder Oil	N-HW	
Speed-Clean by Midmark	N-HW	
Spot Endoscopic Marker	APW	
Stainless Steel Cleaner	HW	Flammable
STARSWAB II	APW	
Static Guard Spray	HW	Flammable
Stavudine	APW	
Sterile Water Injection	APW	
Steris 20 Sterilant Concentrate	HW	Corrosive (Acid)
Storz Endoscope Cleaner	APW	
Streptozotocin	HW	Toxic
Sulfamethoxazole	APW	
Sulfur Precipitated	HW	Flammable
Super Oil Lubricant	N-HW	
Super Sani-Cloth Germicidal Wipes	HW	Flammable
Surepath Liquid Based Pap Test	APW	
SurePath Preservative Fluid	HW	Flammable
Surface Disinfectant & Deoderant	HW	Flammable
Surgical Lubricant Sterile Bacteriostatic	T	
Surgilube Bacteriostatic Lubricant	T	
Swiss Pine	N-HW	
TBS Shur Freeze	N-HW	

Technical Concepts Air Neutralizer	HW	Flammable
Techni-Care Surgical Scrub	D	
Tegaserod Maleate	APW	
Telithromycin	APW	
Telmisartan	APW	
Tempo Ultra SC Insecticide	N-HW	
Ten 20 Conductive EEG Paste	T	
Tenofovir	APW	
Tenofovir Disoproxil Fumarate	APW	
Tenormin	APW	
TexChlor Tablets	N-HW	
TexChlor/Water Solution	D	
Thimerosal preserved medications	HW	Toxic
Thin Prep-Cylolyt	HW	Flammable
Thymol Crystals	T	
Tilex	HW	Corrosive (Base)
Tinidazole	APW	
Tiotropium Bromide	APW	
Tipranavir	APW	
TNT Disinfectant Spray	HW	Flammable
Tolcapone	APW	
Tolterodine Tartrate	APW	
Toluidine Blue O 87% Powder	N-HW	
Topicort Cream	APW	
Topiramate	APW	
Trandolapril	APW	
Triad Adhesive Tape Remover Pads (UNUSED)	APW	
Triad Benzalkonium Chloride Towelettes (UNUSED)	T	
Triad Lubricating Jelly	APW	

Triamcinolone Acetonide	APW	
Trichloroacetic Acid Crystals	N-HW	
Trichloroacetic Acid Solution	HW	Corrosive (Acid)
Trimethoprim	APW	
Tropium Chloride	APW	
Tuberculin Purified Protein Derivative	APW	
Tuberculin Vaccine	APW	
Tums	APW	
Uristix Strips	T	
Urofollitropin	APW	
Vacutainers (UNUSED)	APW	
Valacyclovir	APW	
Valsartan	APW	
Vancomycin	APW	
Vaseline Petroleum Jelly	N-HW	
Vassopressin	APW	
Venlafaxine	APW	
Verapamil	APW	
Verteporfin	APW	
Vesphene II SE - Concentrate	HW	Corrosive (Base)
Vesphene II SE - Diluted Working Solution	N-HW	
Vicriviroc	APW	
Videx & Videx EC	APW	
Vilazodone	APW	
Vinblastine	APW	
Vincristine	APW	
Vitamins - <b>containing Chromium or Selenium</b>	HW	Toxic
Vitamins - <b>not containing Chromium or Selenium</b>	APW	
VSL #3	APW	

Warfarin ≤0.3% tablets	HW	Toxic
Warfarin >0.3% tablets & their packaging	HW	Toxic
WD-40	HW	Flammable
Web Foaming Cleanser with Bleach	N-HW	
White Petrolatum	N-HW	
Windex Liquid	D	
Windex Aerosol	HW	Flammable
Windex Foaming Glass Cleaner	HW	Flammable
Wright Giemsa Stain	HW	Flammable
WYK Safety Sorbent (UNUSED)	T	
Xeroform Gauze	T	
Xeroform Petrolatum Dressing	APW	
Xylocaine	APW	
Zep Glass Cleaner	N-HW	
Zidovudine	APW	
Zileuton	APW	
Zinc Oxide Cream	APW	
Ziprasidone Hydrochloride & Mesylate	APW	
ZNPVA	HW	Flammable/Corrosive (Acid)
Zoledronic Acid	APW	
Zolmitriptan	APW	
Zonisamide	APW	
Spill cleanup debris containing any HW chemicals	HW	
Spill cleanup debris containing any N-HW chemicals	N-HW	

**Key to Abbreviations:**

APW = Anatomical/Pathological Biomedical Waste

D = Drain Disposal

HW = Hazardous Waste

N-HW = Non-Hazardous State Regulated Waste

T = Regular Trash

**Drain Disposal Instructions:**

- Wear chemical resistant gloves, eye protection and a lab coat.
- Flush with at least 10-20 fold excess of water after drain disposal to thoroughly rinse out the sink and to dilute the waste.

[NOTE: For proper disposal of any chemicals not listed above, contact JoAnn Farrell at joann.farrell@yale.edu or 203-432-9384.](mailto:joann.farrell@yale.edu)