

REGISTRATION OF HIGH-SPEED CELL SORTERS

The Yale EHS Biosafety Office has created a [Cell Sorter Registration Form](#) that your lab can utilize to register your cell sorter with EHS. The Biosafety group, along with representatives from the cell sorter biosafety group, will review your application, applicable containment devices that are in place, and your biosafety sorting protocol along with any emergency response standard operating procedures for your sorter. Yale EHS also has sample standard operating procedures and emergency protocols that your lab can use as a starting template to prepare your own “site-specific” biosafety cell sorting documents. Please contact the Yale EHS Biosafety Office or your Yale EHS Safety Advisor to request sample documents.

This group will also collaborate with you to arrange biocontainment testing of your sorter using fluorescent dragon green beads as specified by The International Society for the Advancement of Cytometry (ISAC).